

VOROTEK

experience. knowledge. performance



Portable Suction Pump


EAR WAX REMOVAL BY MICRO SUCTION



Light Weight & Compact



Suction Cannister



Bacterial Filter



Silicone Tubing

Portable Suction Pump

EAR WAX REMOVAL BY MICRO-SUCTION

The VMS-PY Suction Pump is a fully portable electrically powered device that is ideal for effective ear wax removal by micro-suction. It offers a gentle suction flow rate of 15-22L/min and a max vacuum level 75kPa. The VMS-PY Suction Pump is light weight, quiet and fully portable.

Light Weight & Fully Portable

The VMS-PY Suction Pump is fully portable, light weight (3.9kg) with a small compact footprint (314 x 123 x 233mm)

Quiet Operation

The VMS-PY operates quietly at approx. 50 dB sound level.

Optimum Flow Rate

The VMS-PK Suction Pump operates at a maximum flow rate of 15-22 L/min.

Gentle & Effective Vacuum Level

The VMS-PY Suction Pump can be set to operate between 40-75 KPa for effective micro-suction performance.

Simple Operation

The VMS-PY is simple to set up and maintain. Always read the Instructions For Use supplied.

Sturdy Design

The VMS-PY features anti over-flow technology and a shatter - proof & reusable polycarbonate cannister.

Micro-suction Accessories

Vorotek stocks a wide range of accessories required for micro-suction including the VorOtek O Scope for direct vision, suction tips, controllers, tubing & spare parts.

VorOtek Pty Ltd
23 Point Henry Road
Moolap VIC 3224 Australia

T +61 3 5248 5809
E office@vorotek.com.au
www.vorotek.com

VOROTEK
experience. knowledge. performance